

# CUSTOMER SERVICE BULLETIN

## Radio Shack® CUSTOMER SERVICE BULLETIN

400 Atrium—One Tandy Center, Fort Worth, Texas 76102

Bulletin No.	1568-Ø5
Description	Med. Office Systems
Model	Model III
Stock No.	26-1568 Version 1.Ø
Date	Feb. 25, 1982
Page 1 of	2

### Problem

Passwords of less than 8 characters cause files to be opened on drive Ø if the correct disk is not already in place on another drive.

### Solution

The problem is corrected by following the steps listed below.

1. Backup the diskette(s) and make the changes on the Backup copy of the program.
2. In BASIC load the program by typing LOAD"MENU".
3. Make the following corrections:

CHANGES (Retype the line or refer to the Edit section of the owners manual)

Old Line: 59Ø GET1,3:NØ=CVI(NØ\$):DTØ=DAØ:PØ=PZØ:Z1=CVI(Z1Ø):Z2=CVI(Z2Ø):  
Z3=CVI(Z3Ø):DPØ=Z4Ø:BPØ=Z5Ø:YNØ=Z6Ø:CLOSE1:RETURN

New Line: 59Ø GET1,3:NØ=CVI(NØ\$):DTØ=DAØ:PØ=MIDØ(PZØ,1,  
INSTR(PZØ+" "," ")-1):Z1=CVI(Z1Ø):Z2=CVI(Z2Ø):Z3=CVI(Z3Ø):  
DPØ=Z4Ø:BPØ=Z5Ø:YNØ=Z6Ø:CLOSE1:RETURN

Examine lines 135Ø - 136Ø. They should look like this:

```
135Ø DEFFNA(X)=PEEK(X)+PEEK(X+1)*256:E=FNA(&H4ØA4)
136Ø S=FNA(E):IFFNA(S+2)<1ØØØTHENE=S:GOTO136ØELSE
PRINT"TP=TP+"S-FNA(&H4ØA4):END
```

If they do not look like this, add them as shown.

Type 'GOTO135Ø

# CUSTOMER SERVICE BULLETIN

Subject 1568-05 Medical Office Systems Page 2 Of 2

The computer will respond with 'TP=TP+NNNN (a number)

Edit Line 1110 so that the statement TP=TP+38NN equals the number you just got from the computer.

(Line 1110 looks like this before the change.)

```
1110 CLS:PRINT@467,"** LOADING PROGRAM **":GOSUB570:GOSUB580:  
      TP=PEEK(16548)+PEEK(16549)*256:TP=TP+3836:POKE16424,PEEK(16548):  
      POKE16425,PEEK(16549):POKE16549,CINT(TP/256):  
      POKE16548,TP-256*CINT(TP/256):RUNPG$(PR)
```

4. Type SAVE"MENU" to save the changes in the program.
5. At TRSDOS Ready, make a backup copy of the corrected diskette.